

Basic tutorial for FURTHRmind




Create Group



A group contains experiments that share some characteristics, for example the variables of the experiments, and the experiments inside a group differ from each other in respective settings.



Create Group



The screenshot displays the FURTHRmind software interface. At the top, the title bar shows 'FURTHRmind' and 'Projects / Topics Tutorial'. The main interface is divided into several panels: 'Dashboard', 'Comparison', and 'Converge Inspector'. The 'Groups / Samples' panel is active, showing a search bar and radio buttons for 'Groups' (selected) and 'Samples'. A blue button labeled 'Add group' is highlighted by a mouse cursor. The 'Experiments' panel is empty. The 'Diagrams' panel on the right contains a comment field and a 'Raw data' section. The user's email 'daniel.menne@furthr-research.com' is visible in the top right corner.

Right click in the tree view menu in the left and select "Add group"

Create Group



The screenshot shows the FURTHRmind software interface. The main window has a title bar with 'FURTHRmind' and a menu icon. Below the title bar, there's a navigation bar with 'Projects / Topics' and a dropdown menu set to 'Tutorial'. The main area is divided into three panels: 'Groups / Samples', 'Experiments', and 'Diagrams'. The 'Groups / Samples' panel has a search field and radio buttons for 'Groups' (selected) and 'Samples'. The 'Experiments' panel is mostly empty. The 'Diagrams' panel has a search field and a 'Comment:' field. A 'Raw data' section is visible at the bottom right. A dialog box titled 'Add entry - FURTHRmind' is open in the center, with a 'Name:' field containing 'Tutorial' and three radio buttons: 'Blank' (selected), 'import fields from ...', and 'import fields + data from ...'. The dialog has 'Cancel' and 'Ok' buttons at the bottom. A blue arrow points to the 'Ok' button.

Type the group name and
click "Ok"



Create Group



The screenshot shows the FURTHRmind software interface. At the top, the title bar reads 'FURTHRmind' and 'Projects / Topics Tutorial'. The main interface is divided into three panels: 'Groups / Samples', 'Experiments', and 'Diagrams'. In the 'Groups / Samples' panel, there is a search bar, radio buttons for 'Groups' (selected) and 'Samples', and a list containing 'Tutorial'. A blue mouse cursor is pointing at the 'Tutorial' entry. Below the list, there is a 'Name:' field containing 'Tutorial' and an 'Add new field:' dropdown menu. The 'Experiments' panel is a large empty box with a large grey plus sign in the center. The 'Diagrams' panel on the right contains a 'Comment:' text area and a 'Raw data' section with two empty columns.

The first group of this project is created

