

Basic tutorial for FURTHRmind



Create Group by Template



Let's suppose you want to create a group that has the same fields but with different numbers (for example, only the concentration is different). You don't have to repeat the steps above again. Please follow the next slides



Create Group by Template



The screenshot shows the FURTHRmind software interface. The top bar displays 'Projects / Topics' with a dropdown menu set to 'Tutorial'. The user's email 'daniel.menne@furthr-research.com' is visible in the top right. The main interface is divided into three panels: 'Groups / Samples', 'Experiments', and 'Diagrams'. The 'Groups / Samples' panel on the left has a search bar and radio buttons for 'Groups' (selected) and 'Samples'. A dropdown menu shows 'Tutorial' selected, with sub-items 'All experiment' and 'Exp without sample'. Below this, there are input fields for 'Layer Number' (2.0), 'Time' (10.0 s), and 'Concentration' (0.5 mol/l), with the 'Concentration' field highlighted by a blue box. An 'Add new field' section is at the bottom. The central 'Experiments' panel is empty with a large grey plus sign. The right panel contains 'Diagrams' and 'Raw data' sections, both currently empty.

For example, you want to create a group that is the same as “Tutorial”, but only the field “Concentration” is different.



Create Group by Template



The screenshot shows the FURTHRmind software interface. At the top, the title bar reads 'FURTHRmind' and 'Projects / Topics Tutorial'. The main window has three tabs: 'Dashboard', 'Comparison', and 'Convergence Inspector'. The 'Groups / Samples' panel on the left contains a search field, radio buttons for 'Groups' (selected) and 'Samples', and a list of groups. The 'Tutorial' group is selected, and a context menu is open over it with the following options: 'Add group', 'Remove group', 'Add preselection Field', and 'Create sample'. Below the list are input fields for 'Layer Number' (2.0), 'Time' (10.0 s), 'Concentration' (0.5 mol/l), and 'Explanation'. At the bottom of this panel is an 'Add new field' dropdown. The central 'Experiments' area is empty with a large grey plus sign. The right side of the interface has a 'Diagrams' panel with a comment field and a 'Raw data' panel.

As before, right click on the left column and select “Add group”



Create Group by Template



FURTHRmind

Projects / Topics Tutorial

daniel.menne@furthr-research.com

Dashboard Comparison Convergece Inspector

Groups / Samples

Search:

Groups Samples

Tutorial

Name: Tutorial

Layer Number: 2.0

Time: 10.0 s (second)

Concentration: 0.5

Add new field:

Experiments

Diagrams

Comment:

Raw data

Add entry - FURTHRmind

Name: Tutorial 2

Blank import fields from ... import fields + data from ...

Tutorial

Cancel Ok

This time, click on “import fields + data from...”, and then select “Tutorial”



Create Group by Template



The screenshot displays the FURTHRmind software interface. At the top, the title bar shows 'FURTHRmind' and 'Projects / Topics Tutorial'. The main window is divided into several panels: 'Groups / Samples' on the left, 'Experiments' in the center, 'Diagrams' on the right, and 'Raw data' at the bottom right. A modal dialog box titled 'Add entry - FURTHRmind' is open in the center, with the following fields and options:

- Name:** Tutorial 2
- Options:** Blank, import fields from ..., import fields + data from ...
- Content:** Tutorial
- Buttons:** Cancel, Ok

Below the dialog box, the 'Experiments' panel shows a diagram with a central box and lines extending upwards and downwards. The 'Groups / Samples' panel on the left has a search bar with 'Tutorial' entered and a list of items below it. The 'Diagrams' panel on the right has a large empty box with navigation arrows and a 'Comment:' field. The 'Raw data' panel at the bottom right has two empty columns.

Click on "OK"



Create Group by Template



The screenshot shows the FURTHRmind software interface. The top navigation bar includes the text 'Projects / Topics' and a dropdown menu set to 'Tutorial'. The user's email 'daniel.menne@furthr-research.com' is visible in the top right. The main interface is divided into three panels: 'Groups / Samples', 'Experiments', and 'Diagrams'. The 'Groups / Samples' panel on the left contains a search bar, radio buttons for 'Groups' (selected) and 'Samples', and a list of groups including 'Tutorial' and 'Tutorial 2'. Below the list, there are configuration fields for 'Layer Number' (2.0), 'Time' (10.0 s (second)), and 'Concentration' (1 mol/l (mole)), with the 'Concentration' field highlighted by a blue box. An 'Add new field' dropdown is at the bottom. The 'Experiments' panel in the center is empty with a large grey plus sign. The 'Diagrams' panel on the right has a 'Comment' text area and a 'Raw data' section.

A new group with the same fields is created and you can further modified the fields based on your need

